

## CryoArc™ Pallet Transporter (PUR)

The CryoArc Pallet Transporter offering superior thermal performance, it is made from 2.5 to 6 cm polyurethane panels with a density of 48 kg/m<sup>3</sup>. It is optimized to the standard size of a pallet and to the maximum height for air transport.

The CryoArc Pallet Transporter can be delivered to you already assembled and ready for use.

The CryoArc Pallet Transporter can also be delivered disassembled (4 insulated pallets stacked), one person is enough to assemble it:



Position the eutectic plates in the slides provided



Insert the pallet of products directly into the packaging



Close the door, put on the outer cover and you're ready to ship



Permanent stock is available to deal with seasonal peaks



Already assembled and ready for use



Support 3R sustainable program that offer return logistics in place

[sales@cryopak.com](mailto:sales@cryopak.com)

[www.cryopak.com](http://www.cryopak.com)



### United States

551 Raritan Center Parkway  
Edison, NJ 08837

☎ 1.888.827.3393

### Canada


11000 Parkway Blvd  
Anjou, QC Canada H1J 1R6

☎ 1.888.423.7251

### Europe

1 Rue Berthelot  
Maromme, France 76150

☎ +33 2 32 82 59 65



CryoArc™ Pallet Transporter (PUR) have been configured for incredibly user-friendly conditioning and assembly. Both US and Euro pallet options are available.

### Markets & Applications

- Cold chain solutions
- Food & beverage delivery
- Pharmaceuticals & biotechnology products
- Clinical studies/trials
- Any product requiring temperature control

### Standard Sizes

Item #	Product Name	Temperature Range	ID of Payload Box L x W x H (in)	Payload Mass (L)	Duration (hours)
FPU00142	CryoArc Pallet Transporter U.S. and Canada	2°C - 8°C	47 x 39 x 45	1008	120

Item #	Product Name	Temperature Range	ID of Payload Box L x W x H (cm)	Payload Mass (L)	Duration (hours)
FFPU0004	CryoArc Pallet Transporter EU	2°C - 8°C	120 x 100 x 114.56	1008	120

[sales@cryopak.com](mailto:sales@cryopak.com)    [www.cryopak.com](http://www.cryopak.com)



#### United States

551 Raritan Center Parkway  
Edison, NJ 08837

☎ 1.888.827.3393

#### Canada

11000 Parkway Blvd  
Anjou, QC Canada H1J 1R6

☎ 1.888.423.7251

#### Europe

1 Rue Berthelot  
Maromme, France 76150

☎ +33 2 32 82 59 65